

THEORY OF CUSTOM FIELDS in Empower CDS

PART 1 - FIELD TYPES

Field Type

Defines which object the Custom Field refers to and where it appears in Empower

Learn Empower - Save Time for What Truly Matters

Custom Fields

FIELD TYPES

Input Fields

Used to enter additional information or values for further calculations

Field Type

- Sample
- Result
- Peak
- Sample Set
- Component
- Distribution

Custom Fields

FIELD TYPES

Output Fields

Used to perform calculations and evaluations

Field Type

- | | | |
|---|---------------------------------|------------------------------------|
| <input checked="" type="radio"/> Sample | <input type="radio"/> Result | <input type="radio"/> Peak |
| <input type="radio"/> Sample Set | <input type="radio"/> Component | <input type="radio"/> Distribution |

Custom Fields

INPUT FIELD TYPES

Sample

- Applies to a single sample line
- Additional column in a sample set

Examples:
Batch, Storage Conditions

Custom Fields

INPUT FIELD TYPES

Sample Set

- Applies to the entire sample set
- Appears in the Sample Set Information Window

Examples:

Deviation Nr, Labjournal Nr

Custom Fields

INPUT FIELD TYPES

Component

- Applies to each component in the component editor
- Appears in the Component Editor Window as a separate tab next to Value & Purity tabs
- To be used in further calculations

Examples:
Factor, Molecular Weight

Custom Fields

OUTPUT FIELD TYPES

Peak

- Relates to each integrated peak in the chromatogram
- Appears in the Peaks Table by Review

Examples:

Amount/Factor, Evaluation - is the Peak Area > LOD

Custom Fields

OUTPUT FIELD TYPES

Result

- Relates to a single result for the entire chromatogram
- Appears in the Results window by Review

Examples:

Sum of Impurities-Amounts,
Minimum Area

Custom Fields

OUTPUT FIELD TYPES

Distribution

- Relates to each slice of a peak
- GPC/V option only

Empower Questions?

Check out our **trainings & workshops**
or book a **free 30-min consulting session!**



“

Most trainings don't reach the lab.
Our start in the lab.

”

