

SIGNS OF MANUAL INTEGRATION

in Empower CDS

What to look for during data review

PEAK TYPES

In the Peaks Table, each peak has two letters. They describe how the start and end of the peak were integrated.

For example:

- V = Valley
- B = Baseline
- S = Shoulder

Capital letters = automatic integration

Lowercase letters = manual integration

A small detail –

but a very reliable indicator.

2 “MANUAL” COLUMN

You can check the Manual column:

- in the Review Window → Result Window
- or in the Project window (Results tab, via view filter)

If a result was manually integrated, you'll see a checkmark.

Tip: include this column in your reports, so manual changes are always visible.

3 PEAK & RESULT CODES

Some actions are less obvious and can be overlooked:

- **I03** → user manually deleted a peak (Result code)
- **I04** → user manually deleted a drop line (Peak code)

These actions don't always stand out immediately.

The result may be marked as "manual", but the specific change can still go unnoticed.

Tip: You can create Custom Fields for these codes to detect them automatically during review.

Want to improve efficiency of your Empower workflows?

Book a free **25-min Workflow Review**

- Discuss your current processes
- Identify inefficiencies
- Get practical recommendations



German | English

